



Balansposities

NEOX

Balanspositie mogelijk bij zitdiepte i.c.m. FLX2 rugleuning:	20"				22"					24"					26"				
	5"	6"	7"	8"	5"	6"	7"	8"	200x50	5"	6"	7"	8"	200x50	5"	6"	7"	8"	200x50
D40	✓	2,3,4	3,4	4	2,3,4	3,4	4	x	x	3,4	3,4	4	x	x	2,3,4	3,4	x	x	x
D42	✓	✓	2,3,4	3,4	2,3,4	2,3,4	3,4	4	4	3,4	3,4	3,4	4	4	2,3,4	3,4	4	x	x
D44	✓	✓	2,3,4	2,3,4	✓	✓	2,3,4	3,4	3,4	✓	2,3,4	3,4 *	3,4 **	3,4	2,3,4	2,3,4	3,4	4	4
D46	✓	✓	✓	2,3,4	✓	✓	2,3,4	2,3,4	2,3,4	✓	✓	2,3,4	3,4 **	3,4	2,3,4	2,3,4	3,4	3,4	3,4
D48	✓	✓	✓	✓	✓	✓	✓	2,3,4	2,3,4	✓	✓	✓	2,3,4	2,3,4	✓	✓	2,3,4	3,4	3,4
D50	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	2,3,4	2,3,4
D52	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
D54	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
D56	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓

* Balanspositie 2- Medium actief mogelijk vanaf framehoogte 44 cm

** Balanspositie 2- Medium actief mogelijk vanaf framehoogte 47 cm

✓ Alle balansposities mogelijk